

LT-80 PRO

(1,5 LITERS)

Capacity: 1,5 liters
Internal dimensions: 145x132x100 mm
External dimensions: 175x160x260 mm
Weight: 2,5 Kg



ELECTRICAL CHARACTERISTICS

- Power Supply: 220V.
- Ultrasonic power: 80W (160W p-p).
- Heating resistance: 100 W.

ULTRASONIC CHARACTERISTICS

- Ultrasonic power: 80W (160W p-p).
- Working frequency: 40 KHz. Frequency scanning system (sweep system $\pm 2\%$).
- 1 high-performance IBL piezo-electric transducer on steel-titanium.
- Digital control board:
 - Main switch.
 - Ultrasonic time selection: from 0 to 60 min.
 - Temperature selection: from 0° to 80°C.

GENERAL CHARACTERISTICS AND COMPLEMENTS

- Tank made in stainless steel AISI 316.
- External panel made in print-free steel.
- 1 stainless steel basket.
- Stainless steel cover.

GUARANTEE AND QUALITY:

Every equipment manufactured by TIERRA TECH in conformed European Directives about CE certificate. Our machines are due to rules about electric and magnetic compatibility and security of electrical components: Directives: 2004/108/CEE, 2006/95/CEE, EN 55011, EN 60204, EN 61010. This unit is guaranteed for 1 year in case of any manufacture fault.